



MORAN SCIENTIFIC

www.moranscientific.com.au Pty Ltd



- * Professional complete analysis system at cost effective prices.
- * Fast single and multiple full screen imaging. (1 to 50 images)
- * Beam positioning anywhere on image for point analysis or linescan.
- * EDS, WDS and combined X-ray Maps.



- * Fast quantitation with multiple drift correction modes.
- * Linescans with intensity or quant profiles, incredibly fast point analysis and the ability to show or hide elements of interest.
- * Maps collected from Full Spectrum at each pixel, far superior to ROI.
- * View 'Live' images on screen for setting optimum machine conditions.
- * Variable image pixels from 32x32 to 4000x4000 (16Mb images)
- * Image processing and manipulation.
- * Easy setup of palette thresholds for colouring images.
- * Easy image pseudo colouring with fast colour rotation to observe hidden details.
- * Moveable micron bar, magnification and custom labels. (variable fonts and sizes)
- * Simultaneous display of up to 60 images(X-ray, electron or both).
- * Zoomed X-ray elemental images to show concentration distributions.
- * Send results directly to MoranSci Reporter or cut and paste results to other Windows' programs.

and more...